



# THE Comfi™ RANGE



More information about the Comfi™ range of woundcare solutions and helpful information and tools for healthcare professionals can be found on our website at [www.comfi-range.com](http://www.comfi-range.com)

Synergy Health (UK) Ltd, 1 Western Ave, Matrix Park, Chorley, Lancashire, PR7 7NB. UK  
Telephone: +44 (0)1772 299900 Fax: +44 (0)1772 299901 Email: [healthcaresolutions@synergyhealthplc.com](mailto:healthcaresolutions@synergyhealthplc.com)  
[www.comfi-range.com](http://www.comfi-range.com) [www.synergyhealthplc.com](http://www.synergyhealthplc.com)

Comfi, Comfast, MultiStretch, Easywrap, Comfinette, Comfitube, Comfigauz, Comfipression are trademarks of Synergy Health plc. © Synergy Health plc

*As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.*  
SHI5614



The Comfi™ Range of tubular bandages and garments are one of the UK's most popular brands of woundcare solutions.



We pride ourselves in working in partnership with healthcare professionals to deliver great value products, with the assurance that they are manufactured to the highest standards within our UK facility.

Comfi™ has over ten years of experience developing woundcare solutions and our comprehensive range provides solutions for dressing retention, dermatological skin conditions, muscle strains and sprains, and venous disease.

#### Woundcare Solutions

Dressing Retention

Eczema, dry skin and psoriasis

Strains and sprains

Venous Disease

# Comfifast™ Tubular Bandages

Material: Elasticated Viscose

Comfifast™ Tubular Bandages are the preferred product of choice for many clinicians. Versatile, they're commonly used to for dressing retention, wet and dry wrapping of atopic eczema and skin covering for fragile skin. Used for dressing retention, Comfifast™ holds the primary wound dressing in place, providing an extra layer of protection during the healing process.

These colour coded bandages have been produced using light and breathable fabric for patient comfort whilst offering single stretch in the radial direction.

## PRODUCT OVERVIEW

- Holds primary dressing in place without constriction
- No need for pins or tape
- Available in 5 widths and a variety of lengths
- Wash and reuse again and again

## Indicated for...

Dressing Retention

Eczema, dry skin and psoriasis

Covering fragile skin



WIDTH (unstretched)	APPLICATION	LIMB CIRC.	REF 1m	REF 3m	REF 5m	REF 10m
3.5cm Red Line	Small Limbs	8-15cm	F12	-	-	F15
5.0cm Green Line	Medium Limbs	10-25cm	F22	F23	F24	F25
7.5cm Blue Line	Large Limbs	20-45cm	F32	F33	F34	F35
10.75cm Yellow Line	Extra Large Limbs, Head, Trunk (Child)	35-65cm	F42	F43	F44	F45
17.5cm Beige Line	Trunk (Adult)	50-120cm	F52	-	-	F55

# Comfast™ MultiStretch™ Tubular Bandages

Material: Elasticated Viscose

Comfast™ MultiStretch™ Tubular Bandages are our premium range of tubular bandages produced using the same light and breathable fabrics as Comfast™ with added MultiStretch™ technology. Comfast™ MultiStretch™ are able to stretch in the radial, longitudinal and diagonal directions for complete freedom of movement and enhanced patient comfort.

## PRODUCT OVERVIEW

- Holds primary dressing in place without constriction
- No need for pins or tape
- Available in 5 widths and a variety of lengths
- Wash and reuse again and again

## Indicated for...

Dressing Retention

Eczema, dry skin and psoriasis

Covering fragile skin

Primary and secondary lymphoedema

Under compression bandaging



WIDTH (unstretched)	APPLICATION	LIMB CIRC.	REF 1m	REF 3m	REF 5m	REF 10m
3.5cm Red Line	Small Limbs	8-15cm	FMS12	-	-	FMS15
5.0cm Green Line	Medium Limbs	10-25cm	FMS22	FMS23	FMS24	FMS25
7.5cm Blue Line	Large Limbs	20-45cm	FMS32	FMS33	FMS34	FMS35
10.75cm Yellow Line	Extra Large Limbs, Head, Trunk (Child)	35-65cm	FMS42	FMS43	FMS44	FMS45
17.5cm Beige Line	Trunk (Adult)	50-120cm	FMS52	-	-	FMS55

# Comfast™ Easywrap™ Garments

Material: Elasticated Viscose

Comfast™ Easywrap™ Garments are our range of garments for head to toe dressing in the treatment of atopic eczema. The garments are designed to fit each area of the body from a clava for the head through to socks for the feet and provide a skin covering when using steroids, heavy emollients or ointments.

The range has been developed using MultiStretch™ technology and are able to stretch in the radial, longitudinal and diagonal directions for complete freedom of movement and enhanced patient comfort.

The range is completely washable and can be machine washed for frequent use. The garments are manufactured using 'stay soft' fabric to ensure they retain their softness even after repeated washing.

## PRODUCT OVERVIEW

- **Unobtrusive when worn under clothing**
- **Ready to wear, easy to use and quick to apply**
- **Remains in place without the need for additional fixation pins or tape**
- **Wash and reuse again and again**

## Indicated for...

Eczema, dry skin and psoriasis

Covering fragile skin

Dressing Retention

Garments can help in the treatment of atopic eczema;

- **Creates a barrier to scratching**
- **Prevents transfer of creams and emollients to clothing and bedding**
- **Helps support patient concordance**
- **Light fabric allows the skin to breath**



PRODUCT	SIZE: LIMB CIRC.	REF
 Adult Vests	Small - Chest: 71 - 91cm (28" - 36") Medium - Chest: 91 - 117cm (36" - 46") Large - Chest: 117 - 152cm (46" - 60")	CV20 CV21 CV22
 Adult Leggings	Small - Waist / Hip: 71 - 91cm (28" - 36") Medium - Waist / Hip: 91 - 117cm (36" - 46") Large - Waist / Hip: 117 - 152cm (46" - 60")	CL20 CL21 CL22
 Child Vests	6 - 24 months 2 - 5 years 5 - 8 years 8 - 11 years 11 - 14 years	CV24 CV5 CV8 CV11 CV14
 Child Tights	6 - 24 months	CT24
 Child Leggings	2 - 5 years 5 - 8 years 8 - 11 years 11 - 14 years	CL5 CL8 CL11 CL14
 Child Clavas	6 months - 5 years 5 - 14 years	CC5 CC14
 Child Mittens	Up to 24 months (with thumb hole) 8 years (with thumb hole) 8 - 14 years (without thumb hole)	CM24 CM8 CM14
 Child Socks	Up to 8 years 8 - 14 years	CS8 CS14

# Comfinette™, Comfitube™ and Comfigauz™ Tubular Bandages

Material: Comfinette™, 100% viscose; Comfitube™, 100% unbleached cotton; Comfigauz™, 100% bleached cotton

Comfinette™, Comfitube™ and Comfigauz™ are a range of cost effective tubular bandages for dressing retention, designed for use by healthcare professionals. The range is available in 20 metre rolls and 8 sizes suitable to fit all areas of the body from fingers and toes to adult trunks.

**Comfinette™** is made from 100% viscose, ideal for frequent dressing changes.

**Comfitube™** is made from 100% unbleached cotton, containing natural oils making the bandage water repellent.

**Comfigauz™** is made from 100% bleached cotton, ideal for use where natural fibres are preferred.



## PRODUCT OVERVIEW

- **Cost effective**
- **Light and breathable fabrics**
- **Available in a range of sizes to fit all areas of the body**

## Indicated for...

Dressing Retention

SIZE	APPLICATION	COMFINETTE™ 20m	COMFITUBE™ 20m	COMFIGAUZ™ 20m
00	Toes	N006	T006	Z006
01	Fingers and Toes	N016	T016	Z016
12	Bulky Dressings for Fingers and Toes	N126	T126	Z126
34	Children's Limbs / Adult Hands	N346	T346	Z346
56	Adult Limbs	N566	T566	Z566
78	Large Adult Limbs	N786	T786	Z786
T1	Small Trunk / Adult Head	NT16	TT16	ZT16
T2	Large Trunk	NT26	TT26	ZT26

# Comfigrip™ Elasticated Tubular Bandage

Material: Elasticated Viscose

Comfigrip™ Tubular Bandages are our range of elasticated support bandages designed for light support of muscle strains, sprains and joint support. These types of injuries often result in pain, swelling, bruising and loss of movement. Light support to the affected area using Comfigrip™, will support the limb and allow the patient to continue with life as normal while the injury heals.

## PRODUCT OVERVIEW

- Wide range of sizes available to fit all areas of the body
- Provides light support and easy to apply
- Wash and reuse again and again

## Indicated for...

Muscle strains

Sprains

Joint support



WIDTH (unstretched)	APPLICATION	LIMB CIRC.	REF 0.5m	REF 1m	REF 10m
<b>A</b> 4.50cm	Infant Feet and Arm	10-13cm	-	-	GAN5*
<b>B</b> 6.25cm	Hand / Extra Small Ankle	13-16cm	GBN1	GBN2	GBN5*
<b>C</b> 6.75cm	Medium Wrist / Small Ankle and Arm	16-20cm	GCN1	GCN2	GCN5*
<b>D</b> 7.50cm	Large Arm / Medium Ankle	20-24cm	GDN1	GDN2	GDN5*
<b>E</b> 8.75cm	Large Ankle / Medium Knee	24-28cm	GEN1	GEN2	GEN5*
<b>F</b> 10.00cm	Large Knee / Medium Thigh	28-36cm	GFN1	GFN2	GFN5*
<b>G</b> 12.00cm	Large Thigh	36-46cm	GGN1	GGN2	GGN5*

\*Comfigrip™ Latex Free option is also available for patients allergic to natural rubber latex.

# Comfipression™ Hosiery Kits

Material: Elasticated Viscose

Comfipression™ Hosiery Kits are designed to treat and manage venous disease to provide compression therapy to the affected limb. Comfipression™ is ideal for patients whose healthcare professional has given them the option of wearing compression hosiery instead of traditional compression bandages.

**Hosiery Kit;** 2-piece compression kit containing;  
2 x first layer compression liners 10mmHg  
1 x second layer compression stocking 30mmHg

**Liner Kit;** replacement liner kit containing;  
2 x first layer compression liners 10mmHg.



## PRODUCT OVERVIEW

- Discreet and comfortable to wear under clothing and with shoes
- Consistently applies a known pressure of 40mmHg compression at the ankle
- Safe and easy to use
- Wash and reuse again and again up to 100 washes

## Indicated for...

Venous disease

Leg ulcers

Varicose veins

Managing oedema

Controlled lymphoedema

SIZE	THREAD COLOUR	HOSIERY KIT	LINER KIT
Small	Red Line	CPHK01	CPLK01
Medium	Blue Line	CPHK02	CPLK02
Large	Green Line	CPHK03	CPLK03
X Large	Yellow Line	CPHK04	CPLK04
XX Large	Purple Line	CPHK05	CPLK05