



www.koomedical.com



KOO MEDICAL EQUIPMENT

KOO MEDICAL EQUIPMENT

Koo Medical Equipment is a multinational company that has been working for many years in the medical industry. Its headquarters and its production unit are located in the Far East.

The company specialises in the manufacture of plastics for medical use.

Our staff's commitment and expertise are a guarantee in the development and production of specialized products used for anesthesia or respiratory support.

Koo Medical Equipment has branches in Europe and Asia.

Every branch provides sales and, marketing support and constant group's collaboration allows the development of innovative solutions in the design and manufacture of medical devices.

A continuous control of raw materials used in the manufacture of our products and of the services offered to our partners is the basis of all the activities carried out by Koo Medical Equipment.

Each product designed and manufactured by Koo satisfies the strictest qualitative standards required by the market and all products are certified according to the European standards for Medical Devices attested by DNV-BUSINESS ASSURANCE-NORWAY-CE0434 class IIa.

The production plants are certified ISO 9001 and ISO 13485 and FDA approval.

The respect for absolute quality in the anesthesia mask's production makes the company a world class manufacturer.

The same attention is used everyday in all the activities carried out by the company in order to maintain high-quality standards both for the devices manufactured and for the services offered, thus achieving an excellent relationship with our partners.

Every product is manufactured with PVC materials, DEHP and BPA free

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



DNV BUSINESS ASSURANCE
EC CERTIFICATE - FULL QUALITY ASSURANCE SYSTEM

Certificate No. 75900-2010-CE-RGC-NA 1.0
 This Certificate consists of 3 pages

This is to certify that the Quality Management System of

Koo (Shanghai) Industries Co., Ltd.

100 Zhongde Rd, Songjiang Shanghai 201614, P. R. China

for design, production and final product inspection/testing of
**Sterile and Non-sterile Medical Devices for Anesthesia, Respiratory
 and Critical Care**

has been assessed with respect to
 the conformity assessment procedure described in Article 11.3.a and Annex II excluding section 4
 (Module H) of Council Directive 93/42/EEC on Medical Devices, as amended, and found to comply

Further details are given overleaf

Place and date:
 Hovik, 19 November 2013

For DET NORSKE VERITAS CERTIFICATION AS
 NORWAY

Aud Loken Eiklid

Aud Loken Eiklid
 Certification Manager



Notified Body No.:
 0434

This Certificate is valid until:
 25 May 2015

Alexey Shiryaev
 Technical Reviewer

This Certificate has been digitally signed. See www.dnv.com/digital/signatures for more info

Notice: The certificate is subject to terms and conditions overleaf. Any significant changes in design or construction may render this certificate invalid.
 If any product, system or service is found to be non-compliant, the certificate holder shall be notified as soon as possible. The certificate holder shall be notified as soon as possible. The certificate holder shall be notified as soon as possible.



DNV BUSINESS ASSURANCE
MANAGEMENT SYSTEM CERTIFICATE

Certificate No. 47446-2009-AQ-RGC-NA

This is to certify that the Management System of:

Koo Medical equipment (Shanghai) Co., Ltd.

1 Songjiang Industrial Park, Songjiang, Shanghai, 201613 P. R. China

has been found to conform to the standard:

ISO 13485:2003

This Certificate is valid for the following product or service ranges:

**Design and manufacture of both Sterile and Non-sterile medical devices for
 Anesthesia, Respiratory and Critical Care.**

Place and date:
 Hovik, 28 February 2012
 for the Accredited Unit:
 DET NORSKE VERITAS
 CERTIFICATION AS, NORWAY

Eugenie Winger Husebye

Eugenie Winger Husebye
 Management Representative



As set out in the Certification Agreement may render this Certificate invalid.
 If this has been digitally signed. See www.dnv.com/digital/signatures for more info.
 DET NORSKE VERITAS CERTIFICATION AS, VERITASVEIEN 1, 1322 HOVIK, NORWAY. TEL: +47 67 57 99 00 FAX: +47 67 57 99 11 - www.dnv.com



DNV BUSINESS ASSURANCE
MANAGEMENT SYSTEM CERTIFICATE

Certificate No. 47453-2009-AQ-RGC-NA

This is to certify that the Management System of:

Koo Medical equipment (Shanghai) Co., Ltd.

1 Songjiang Industrial Park, Songjiang, Shanghai, 201613 P. R. China

has been found to conform to the standard:

ISO 9001:2008

This Certificate is valid for the following product or service ranges:

**Design and manufacture of both Sterile and Non-sterile medical devices for
 Anesthesia, Respiratory and Critical Care.**

Initial Certification date:
 20 December 2000

This Certificate is valid until:
 7 March 2015

The audit has been performed under the
 supervision of

Jin Ming Xu
 Lead Auditor



Place and date:
 Hovik, 28 February 2012
 for the Accredited Unit:
 DET NORSKE VERITAS
 CERTIFICATION AS, NORWAY

Eugenie Winger Husebye

Eugenie Winger Husebye
 Management Representative

Lack of fulfillment of conditions as set out in the Certification Agreement may render this Certificate invalid.
 This Certificate has been digitally signed. See www.dnv.com/digital/signatures for more info.
 HEAD OFFICE: Det Norske Veritas AS, Veritasveien 1, 1322 Hovik, Norway. Tel: +47 67 57 99 00 Fax: +47 67 57 99 11 - www.dnv.com



INDEX

Page

Anaesthesia Masks	4-5
Laryngeal Masks	6-7
CPAP Masks and Peep Valve	8-9
Non-invasive ventilation Masks NIV	10-15
Accessories for CPAP and NIV Masks	16
Manual resuscitators	17-18
Respiratory exercise systems	19
Oxygen masks	20-22
Nasal cannula	22
Connection tubes for oxygen therapy and aerosol therapy	23
Aerosol masks and nebulizer kits	24-25
Fisiojet nebulizer/compressor system	26
Fisiojet nebulizers	26-27
Fisiojet masks and accessories for aerosol therapy	28
FISIOCHAMBER VISION Spacer	29-30
Spacer masks	30
Veterinary - masks and accessories	31-32

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



www.koomedical.com

Anaesthesia mask

Anaesthesia Mask CRYSTAL™



- Inflating air cushion with inflation valve.
- It fits the patient's face and optimizes the application of anesthetics.
- PVC - DEHP and BFA free high transparency allows good visibility.
- Latex-free
- Safe and easy to use.
- Single patient use.
- Available in 8 sizes with color-coded fixing rings.

code	description	size	weight gr	dead space	colors	connector	items per box
KM-200	anaesthesia mask	extra large	59,5	310 ml	transparent	22F	50
KM-201	anaesthesia mask	large	55,5	250 ml	red	22F	50
KM-202	anaesthesia mask	medium	51,0	220 ml	yellow	22F	50
KM-203	anaesthesia mask	small	37,5	180 ml	green	22F	50
KM-204	anaesthesia mask	pediatric	27,0	85 ml	white	22F	50
KM-205	anaesthesia mask	infant	17,5	45 ml		15M	50
KM-206	anaesthesia mask	neonatal	14,5	25 ml		15M	50
KM-207	anaesthesia mask	pre-neonatal	11,0	20 ml		15M	50

- Each mask is packed individually with a description of the contents.
- Appreciated in each country; they guarantee safety and comfort.
- Vanilla scent.

by request they can be delivered with a strawberry scent with the same code and adding the letter S= KM-200S, KM-201S, KM-202S, KM-203S, KM-204S, KM-205S, KM-260S, KM-207S



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Anaesthesia Mask

CRYSTAL™ PLUS

- Comfortable inflating cushion with inflation valve.
- PVC - DEHP and BFA free high transparency allows good visibility.
- Latex-free
- Safe and easy to use.
- Patented anaesthesia mask designed with lateral grips for a comfortable and confident mask manipulation.
- Single patient use.
- Available in 7 sizes with fixing rings.



code	description	size	weight gr	dead space	connector	items per box
KM-279	anaesthesia mask	large	55,5	250 ml	22F	50
KM-278	anaesthesia mask	medium	51,0	220 ml	22F	50
KM-277	anaesthesia mask	small	37,5	180 ml	22F	50
KM-276	anaesthesia mask	pediatric	27,0	85 ml	22F	50
KM-275	anaesthesia mask	infant	17,5	45 ml	15M	50
KM-274	anaesthesia mask	neonatal	14,5	25 ml	15M	50
KM-273	anaesthesia mask	pre-neonatal	11,0	20 ml	15M	50

- Each mask is packed individually with a description of the contents.



“CRYSTAL PLUS is considered to be the best mask for anaesthesia. Comfortable and confident manipulation.”

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



www.koomedical.com

Laryngeal mask

Laryngeal mask

Crystal Airway Mask™



- In PVC -DEHP and BFA free.
- Latex-free.
- Soft and non-traumatic laryngeal protector, avoiding lesions and ensuring a good hold and ventilation.
- The rear cuff is reinforced with a thicker wall so that the device can be easily inserted.
- Inflation line with valve is integrated in to the wall of the ventilation tube.
- Sterile.
- Single patient use.
- Available in 7 sizes. Each size is color coded.

code	size	patient's weight (kg)	diameter (mm) I.D. internal	diameter (mm) O.D. external	maximum volume of the cuff (ml)	items per box
KM-877	1	< under 5	5,6	8,6	5	10
KM-867	1,5	5-10	6,4	10	7	10
KM-857	2	10-20	7,8	12	10	10
KM-847	2,5	20-30	9,5	14	15	10
KM-837	3	30-50 small adults and children	10	16	20	10
KM-827	4	50-70 men-women	11,4	18	35	10
KM-817	5	> over 70 adults/men	11,9	20	48	10



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com


Laryngeal mask

FLEXCAM reinforced

- In PVC - DEHP and BPA free with a metallic reinforcement.
- Latex-free.
- Soft and non-traumatic cuff, avoiding lesions and ensuring a good hold and ventilation.
- The rear cuff is reinforced with a thicker wall so that the device can be easily inserted.
- Inflation line with valve is integrated in to the wall of the ventilation tube.
- Sterile.
- Single patient use.
- Available in 5 sizes. Each size is color coded.



code	size	patient's weight (kg)	items per box
KM-858	2	10-20	10
KM-848	2,5	20-30	10
KM-838	3	30-50	10
KM-828	4	50-70	10
KM-818	5	> over 70	10



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



Ventilation Mask CPAP Zephir™



- Safety valve.
- Inflating cushion with inflation valve.
- It fits perfectly ensuring a high hold and avoiding air leak.
- Provided with a port for oxygen and monitoring.
- In PVC -DEHP and BPA free high transparency allows good visibility.
- Latex-free.
- Safe and easy to use.
- Provided with two ports: Inhaling (higher) port 22M and exhaling (lower) port (36 mm OD x 33 mm ID).
- Single patient use.
- Available in 3 sizes.

code	description	size	weight (gr)	dead space	items per box
KM-301A	mask for CPAP	large	67,5	245 ml	22
KM-301DA	mask for CPAP	large deep	74,5	275 ml	22
KM-303A	mask for CPAP	small	68,5	215 ml	24

- Each mask is packed individually with a description of the contents.



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Valve Peep

- Adjustable from 5 to 20 cm H2O for CPAP ventilation.
- With connection OD 33 mm for CPAP masks KM-301A, KM-301DA, KM-303A.
- With two reducer connectors: 1x22Mx30,5 mm ; 1x22Mx33mm they allow a universal use of the valve for helmets, masks and circuits.
- Single patient use.



- Code: KM-809 (50 items per box).
- Each valve is packed individually with a description of the contents.

Headgear for NIV and CPAP masks



- In elastic and breathable material. 100% acrylic
- Soft and easy to use
- Single patient use

code	description	size	items per box
KM-P007	fixing mask for mask CPAP with two ports	large	100
KM-P008	fixing mask for mask CPAP with two ports	small	100

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



Non-Invasive Ventilation Face mask

Bluestar™

For lung ventilation devices. Modality CPAP and BILEVEL assisted and/or supervised.



- Inflating cushion adaptable to the patient's face with a firm and comfortable fit. No skin breakdown.

The kit is provided with:

PVC-DEHP-BPA free and transparent mask, with inflatable cushion, valve with 6 grips to fix the mask, adapter connector (KM-1019) with ports for O₂ and CO₂, pressure monitoring and a double caps for the introduction of a nasal cannula.

- Pump to inflate and deflate the cushion.
- Headgear's cuff with a rapid removal system.
- Latex Free
- Possible link with connector (KM-903 not included) provided with a hole for the CO₂ discharge and a port for the improvement of O₂.
- Single patient use.

The kit is inside a box with a label that indicates the contents.

code	description	size	weight (gr)	dead space	diameter of the connector	items per box
KPM-001L	face mask NIV	extra large	91,0	290 ml	22F	50
KPM-001M	face mask NIV	large	79,5	240 ml	22F	50
KPM-002M	face mask NIV	medium	70,0	220 ml	22F	50
KPM-003M	face mask NIV	small	65,0	175 ml	22F	50
KPM-004M	face mask NIV	pediatric	46,0	85 ml	22F	50
KM 1019	Adapter connector	standard	19,0	17 ml	22M	100



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Non-Invasive Ventilation nasal mask

Moonlight™

for CPAP BILEVEL non-invasive lung ventilation devices.

- Comfortable and easy to place
- In transparent PVC- DEHP and BFA free, with soft and adaptable flange
- Latex free
- With port and cap for enrichment with O₂ and monitoring
- Rotating connector 360°
- Spacer for correct placement
- Cushions for front support
- Single patient use



code	description	size	weight (gr)	dead space	diameter of the connector	items per box
KM-501	nasal mask NIV	large	53,8	90 ml	22mmF	50
KM-502	nasal mask NIV	medium	52,8	85 ml	22mmF	50
KM-503	nasal mask NIV	small	47,6	75 ml	22mmF	50

Every mask is packed individually with a description of the contents.

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



www.koomedical.com

NIV masks

Non-Invasive Ventilation Face mask **Bluestar PLUS™**

for pulmonary ventilation's devices, modality CPAP and BILEVEL assisted and supervised.



- Innovative design that guarantees an optimal air-leak control and a high patient comfort level.
- Lightweight.
- Soft replaceable cushion minimizes skin breakdown.
- The kit has all the accessories for positioning.

kit contents:

- High transparent PVC- DEPH and BPA free mask with inflating and replaceable cushion, with inflation valve.
- 8 grips for fixing the headgear.
- Frontal adjustable cushion with 2 grips.
- Ports for O₂ and CO₂, pressure's monitoring and a port for introducing gastric tubing and with the possibility of introducing a bronchoscope.
- Two connectors 22M and 22F.
- Pump for inflating and deflating the cushion.
- Headgear with a rapid removal system.
- Latex free
- The mask's connector can be connected to expiratory system (KM-903 not in the kit) provided with a hole for the CO₂ discharge and with a port for the enrichment of O₂.
- Every kit is delivered in a box with a label indicating the contents.



"the best handy and low-cost solution for non-invasive ventilation with a device that allows the replacement of the cushion and of the shell"

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



Non-invasive ventilation face mask

Bluestar PLUS™

code	description	size	weight (gr)	dead space	diameter of the connector	kits per box
KM-310B	face mask NIV	extra large	158	205 ml	22F/M	50
KM-311B	face mask NIV	large	150	195 ml	22F/M	50
KM-312B	face mask NIV	medium	140	165 ml	22F/M	50
KM-313B	face mask NIV	small	140	155 ml	22F/M	50
KM-314B	face mask NIV	pediatric	140	95 ml	22F/M	50

Item KM-310B is provided with mask shell Large and cushion Extra Large.

Inflatable and replaceable cushions for

Bluestar PLUS™

code	description	size	for mask:	items per box
KM-3101	cushion replacement	extra large	KM-310	50
KM-3111	cushion replacement	large	KM-311	50
KM-3121	cushion replacement	medium	KM-312	50
KM-3131	cushion replacement	small	KM-313	50
KM-3141	cushion replacement	pediatric	KM-314	50

Every cushion is packed individually. Mask shells available.



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



www.koomedical.com

NIV masks

Non-invasive ventilation nasal mask

Deluxe™

for of CPAP BILEVEL non-invasive lung ventilation devices



- High transparency PVC - DEHP and BPA free
- Latex-free
- Replaceable cushion with valve
- Comfortable positioned optimally and eliminates air leakage
- New design guarantees even lower skin breakdown
- Ports with caps for enrichment with O₂ and monitoring
- Comfortable frontal cushion
- 360 rotating connector
- Single patient use

Delivered in a kit with positioning accessories.

kit contents

- Mask with cushion, inflating pump, headgear fixing-mask

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

NIV masks



www.koomedical.com



Non invasive ventilation nasal masks and replacement cushion **Deluxe™**

for CPAP BILEVEL non-invasive lung ventilation devices

code	description	size	weight (gr)	dead space	internal diameter's connector	kits per box
with no holes for CO₂ exhalation:						
KM-510	nasal mask NIV	large	71	80 ml	22mmF	50
KM-511	nasal mask NIV	medium	71	75 ml	22mmF	50
KM-512	nasal mask NIV	small	71	70 ml	22mmF	50
with holes for CO₂ exhalation						
KM-513	nasal mask NIV	large	71	80 ml	22mmF	50
KM-514	nasal mask NIV	medium	71	75 ml	22mmF	50
KM-515	nasal mask NIV	small	71	70 ml	22mmF	50



code	description	size	for mask:	items per box
KM-5101	Cushion replacement	large	KM-510-KM-513	50
KM-5111	Cushion replacement	medium	KM-511-KM-514	50
KM-5121	Cushion replacement	small	KM-512-KM-515	50

Every mask is packed individually with a description of the contents.

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



www.koomedical.com

Accessories for CPAP and NIV masks

Headgear fixing-mask for NIV and CPAP masks

Accessory



- In elastic and breathable material. 100% acrylic
- Soft and easy to use
- Single patient use



code	description	size	items per box
KPM-P001	fixing mask NIV	large	100
KPM-P002	fixing mask NIV	small	100
KM-P003	fixing mask NIV	large	100
KM-P004	fixing mask NIV	small	100
KM-P007	fixing mask for mask CPAP with two ports	large	100
KM-P008	fixing mask for mask CPAP with two ports	small	100



NIV masks' connector with a hole for the CO₂ emission and a port for the enrichment of O₂ or monitoring.
Diameter O.D. 22 mm M- I.D. 15 mm F

code	items per box
KM-903	50

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Manual Resuscitator bag

- PVC and PC, DEHP and BFA free resuscitator bag
- Latex free
- For adults 1700ml volume
- Pediatric 700ml volume
- Infant 300ml volume
- Single patient use



kit in a case containing: resuscitator bag, oxygen tube, reservoir, and Crystal mask

code	description	size	with mask:	items per box
KM-806	manual resuscitator	adult	KM-202 medium	10
KM-807	manual resuscitator	pediatric	KM-204 pediatric	10
KM-808	manual resuscitator	infant	KM-205 infant	10



kit in a case containing: resuscitator bag, oxygen tube, reservoir, and Crystal mask, with overpressure pop off control valve

code	description	size	with mask:	items per box
KM-806A	manual resuscitator	adult	KM-202 medium	10
KM-807A	manual resuscitator	pediatric	KM-204 pediatric	10
KM-808A	manual resuscitator	infant	KM-205 infant	10

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



www.koomedical.com

Manual resuscitator

Kit in a case containing: resuscitator bag, oxygen tube, and reservoir, without mask, with overpressure "pop off" control valve

code	description	size	items per box
KM-8062	manual resuscitator	adult	10
KM-8072	manual resuscitator	pediatric	10
KM-8082	manual resuscitator	infant	10



Kit in a case containing: resuscitator bag, oxygen tube, reservoir, without mask, with overpressure "pop off" control valve, and PEEP valve adapter

code	description	size	items per box
KM-8061	manual resuscitator	adult	10
KM-8071	manual resuscitator	pediatric	10
KM-8081	manual resuscitator	infant	10



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Tri-Gym

Lung exercises

- In polycarbonate, balls in PVC DEHP and BFA free, tube and mouthpiece in polyethylene
- Latex free
- For inspiration and expiration
- Two adjustable valves controlling independently inspiration and expiration
- With three measuring cylinders for the control of the flow
- Flow control is done through the three colored balls, which indicate the amount of flow per second 600 cc green, blue 900 cc, 1200 cc red.
- With corrugated tube and mouthpiece
- Single patient use



Code KM-805

a single packaging, 24 items per box.

PipeP kit

System for lung exercises during inspiration and expiration

- In polycarbonate, connectors and resistors in ABS
 - Latex free DEHP and BPA free
- Kit containing:
- One two way valve:
IN= inspiration OD 22mm
OUT= expiration OD 22mm ID 15mm
 - 8 resistors for PEP or PIP with color code with the following holes' diameters: 1,5 mm black, 2,0 mm white, 2,5 mm yellow, 3,0 mm blue, 3,5 mm green, 4,0 mm orange, 4,5 mm pink, 5,0 mm brown; with connector for manometer, OD 15 mm
 - Mask not included in the kit
 - Single patient use



Code KM-803

Every kit is in a box, 24 kits in a carton.

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



www.koomedical.com

Oxygen therapy

Evergreen oxygen masks



- Transparent PVC, DEHP and BFA free
- Latex free
- Very light
- With wide elastic band and adjustable clip
- Single patient use
- Single packaging

code	size	items per box
KM-210	adult	50
KM-211	pediatric	50
KM-214	infant	50

- With 210 cm oxygen supply tubing, kink resistant, with a conical standard connector
- Single packaging
- Single patient use

code	size	items per box
KM-212	adult	50
KM-213	pediatric	50
KM-2141	infant	50



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Evergreen High concentration oxygen masks

- In high transparent PVC - DEHP and BFA free
- Latex-free
- With wide elastic band and adjustable nose clips
- With reservoir and one-way valve
- With 210 cm oxygen supply tubing, kink resistant with a conical standard connector
- Single packaging
- Single patient use



code	size	items per box
KM-260	adult	50
KM-270	pediatric	50

Venturi Evergreen oxygen masks

- In high transparent PVC- DEHP and BPA free
- Latex-free
- With wide elastic band and adjustable nose clips
- VENTURI system with regulator to control the flow of oxygen (24%,26%,28%,31%,35%,40%,50%.)
- With corrugated tube 15 cm
- With 210 cm oxygen supply tubing, kink resistant, with a conical standard connector
- Single packaging
- Single patient use



code	size	items per box
KM-131	adult	50
KM-133	pediatric	50

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



www.koomedical.com

Oxygen therapy

Evergreen tracheostomy oxygen masks



- In high transparent PVC DEHP and BPA free
- Latex free
- Fixing elastic band
- Single packaging
- Single patient use

code	size	items per box
KM-908	adult	50
KM-918	pediatric	50



Evergreen oxygen nasal cannula

- In high transparent PVC DEHP and BPA free
- Latex free
- With 210 cm oxygen supply tubing, with a conical standard connector
- Single packaging
- Single patient use

code	size	items per box
with soft and round nasal tip		
KM-218	adult	50
KM-222	pediatric	50
KM-224	infant	50
KM-226	neonatal	50
with soft and flared tip		
KM-228	adult	50
KM-230	pediatric	50



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Tubes for oxygen and aerosol therapy

- In high transparent PVC - DEHP and BPA free
- Latex-free
- Single packaging

code	size	items per box
with standard conical connector		
KM-2101	210 cm	50
KM-2111	420 cm	50
KM-2106	750 cm	25
with universal connector		
KM-2102	210 cm	50
KM-2103	420 cm	50



Universal Connector



Standard Conical Connector



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



www.koomedical.com

Aerosol therapy

Evergreen Aerosol Masks



- In high transparent PVC- DEHP and BPA free
- Latex free
- Ultra light
- With wide elastic band and adjustable nose clips
- Single packaging
- Single patient use

code	size	items per box
KM-220	adult	50
KM-221	pediatric	50

Nebulizer kit for aerosol therapy



- Evergreen mask in high transparent PVC- DEHP and BPA free
- Latex free
- With wide elastic band and adjustable nose clips
- Nebulizer with certification, tubing 210 cm
- Single packaging
- Single patient use

code	size	items per box
with tubing, standard conical connector		
KM-113A	adult	50
KM-114A	pediatric	50
with tubing, universal connector		
KM-113	adult	50
KM-114	pediatric	50



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Nebulizer kit for aerosol therapy

- With a corrugated tube with a mouthpiece
- Nebulizer with certification
- Tubing, 210 cm
- In high transparent PVC - DEHP and BPA free
- Latex free
- Single packaging
- Single patient use



code	size	items per box
with tubing, standard conical connector		
KM-111	One size	50
with tubing, universal connector		
KM-112	One size	50



Nebulizer high performance for aerosol therapy

- In PVC - DEHP and BPA free
- Latex free
- Single packaging

code	size	items per box
KM-117	One size	50

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



www.koomedical.com

Aerosol **FISIOJET** kit with:



- High quality compressor
- Nebulizer Fisiojet KM 1000 with high performance in polycarbonate
- Latex free DHEP and BPA free
- Elbow connector
- Mouthpiece
- Nose piece
- Connection tube
- In a carry bag

Every item is in a box

Code:

FJ001 Junior with Fisiojet infant Mask and Fisiojet pediatric Mask

FJ002 Senior with Fisiojet Adult Mask



FISIOJET nebulizer for aerosol therapy

- High efficiency nebulizer with inspiratory and expiratory valves
- This type of supply optimizes the inhalation therapy, guaranteeing, maximum medication absorption during inspiration and minimum medication loss during expiration
- Minimal flow 3,0 l/min maximum flow 6,0 l/m
- Minimum cap. 2 ml maximum cap. 8 ml
- In polycarbonate DHEP and BPA free
- Latex free

Provided with the following items:
elbow connector, mouthpiece, nose piece

- single patient use

Code:

KM-1002 in a box



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com
Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

FISIOJET JUNIOR

nebulizer for aerosol therapy with:

- Nebulizer, elbow connector, mouthpiece, nose piece
- Fisiojet infant mask, Fisiojet pediatric mask
- Connecting tube ca. 120 cm
- Single patient use

Code:

KM-1003 in a box



FISIOJET SENIOR

nebulizer for aerosol therapy with:

- Nebulizer, elbow connector, mouthpiece, nose piece
- Fisiojet adult mask
- Connecting tube ca. 120 cm
- Single patient use

Code:

KM-1004 in a box



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



FISIOJET

facial mask for aerosol therapy



- In PVC DEHP and BPA free
- Latex free
- Soft and comfortable
- Minimum drug loss
- A valve for the exhalation of CO₂
- Single packaging
- Single patient use

code	size	items per box
KM-1005	adult	50
KM-1006	small pediatric	50
KM-1007	large pediatric	50
KM-1008	infant	50

Accessories:

code	items per box
KM-1009 Fisiojet nebulizer	50
KM-1010 standard elbow connector	50
KM-1011 nose piece to use with the mouthpiece	50
KM-1012 standard mouthpiece	50
Km-1013 tube of aerosol ca. 120 cm	100

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for EMEA (Europa Middle East Africa) Markets - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for Asian Pacific Market - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

FISIOCHAMBER VISION SPACER

*Designed to be the best therapy
Anti-static Inhaler chamber
with Vision valve for metered
dose Inhaler
BPA free*

- Anti static Chamber: the special antistatic treatment of the chamber ensures a complete suspension of the drug, improving the dose available for the patient.
- Vision valve: the movement of the Vision Valve indicates inhalation and permits monitoring of patient compliance. The one-way valve has very little resistance to flow and can be opened with minimal inspiratory effort (suitable for patients with low tidal volume).
- Universal adapter PROPER FIT: the soft and silicon material of the adapter ensures a perfect seal with all the inhaler connector available on the market.
- Spiral vents: the spiral vent is designed to direct the drug to the one way valve, thereby increasing the quantity of respirable particles
- Mask COMFORT & SEAL: the soft material and the ergonomic contour ensures a perfect seal and comfort to the patient in order to allow for a correct inhalation technique. The safety expiratory valve ensures a complete CO₂ exhalation.
- Fio Chamber Vision is produced with materials free of BPA in compliance with the new regulations.

Materials:

Antistatic inhaling chamber in polycarbonate
Silicone ring valve for universal connection
ABS cap
Safety silicone breathing valve
PVC mask
Latex free DEHP and BPA free



*Antistatic spacer chamber for inhaling
FISIOCHAMBER with VISION valve for measured inhaling spray
single packaging, 12 items per box*

Code: KM-1020

KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com



www.koomedical.com

FISIOCHAMBER vision spacer

FISIOCHAMBER VISION spacer



Antistatic spacer chamber for inhaling
FISIOCHAMBER with VISION valve for measured inhaling spray

single packaging, 10 items per box

Code: KM-1021 - with pediatric mask COMFORT & SEAL, for children aged 3-6

Antistatic spacer chamber for inhaling
FISIOCHAMBER with VISION valve
for measured inhaling spray

single packaging, 10 items per box

Code: KM-1022 - with infant mask
COMFORT & SEAL, for children aged 0-3



FISIO COMFORT & SEAL, mask



*in soft material, ergonomic profile, with safety exhaling
valve that guarantees a complete CO₂ exhalation,
adjustable to FISIOCHAMBER and aerosol therapy items.*

single packaging, 50 items per box

Code: KM-1023 - for adults

Code: KM-1024 - pediatric for children aged 3-6

Code: KM-1025 - infant for children aged 0-3



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com
Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Veterinary mask



- Mask for anesthesia, ventilation and aerosol therapy for dogs, cats and small animals.
- Equipped with a black silicone diaphragm and a highly transparent PVC cone.
- Safe and easy to use.
- Latex-free
- Single packaging



6 sizes available

code	size	weight (gr)	connector		items per box
			internal diameter	external diameter	
KM-701	6 for big dog	122	27 mm	22 mm	25
KM-702	5 for medium dog	115	27 mm	22 mm	25
KM-703	4 for small dog	47	27 mm	22 mm	25
KM-704	3 for big cat	23,5	27 mm	22 mm	25
KM-705	2 for medium cat	13,5	15 mm	12 mm	25
KM-706	1 for small cat	10	15 mm	12 mm	25

The sizes 3,4,5,6 are delivered with a fixing mask ring



www.koomedical.com

Veterinary mask

Diaphragm replacement black silicone for veterinary mask

single packaging

6 sizes available



code	size	items per box
KM-711	for mask size gr. 6	10
KM-712	for mask size gr. 5	10
KM-713	for mask size gr. 4	10
KM-714	for mask size gr. 3	10
KM-715	for mask size gr. 2	10
KM-716	for mask size gr. 1	10

Fixing mask for veterinary, in transparent PVC
single packaging

Code: KM-721 - 10 items per box



KOO MEDICAL EQUIPMENT

Web: www.koomedical.com

Responsible for **EMEA (Europa Middle East Africa) Markets** - Koo Europe s.r.l. viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market** - Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Responsible for **EMEA (Europa Middle East Africa) Markets**
Koo Europe s.r.l.

viale delle Industrie 5 - 20020 Arese (Milano) ITALY Tel: 0039 02 93589495 - Email: info@kooeurope.com

Responsible for **Asian Pacific Market**

Koo Asia Ltd Unit 23,9 Floor Goodwill Ind. Building, No.36-44, Pak Tin Pa Street, Tsuen Wan, Hong Kong
Tel:+852 2110 4780 Fax:+852 2110 4783 - Email: sales@koo-asia.com

Koo Medical Equipment (Shanghai) Co.,Ltd.

www.koomedical.com

1 Songjiang Industrial Park - Songjiang, Shanghai, 201613 P.R.China

