



## Parametric Check pT-Tester

**PATENTED**

The first, double check, indicator for two, most critical parameters on your sealing device, namely temperature and pressure.




pT-Tester is an unique dual parameter check for routine and daily monitoring of your sealer device performance, to ensure complete integrity of your sterile barrier system (SBS).

pT-Tester complies to the ISO EN 11607-2 requirements.

ARROWPACK uses only high developed nano-ink technologies, which allows precise interpretation of parameters.


pT-Tester is usable in conjunction with ARROWPACK Push DyeTest and will ensure complete integrity of all seal areas.

pT-Tester allows CSSD managers to check at any time and without any doubt the pressure and temperature of the sealing devices in a quick and professional way.



## Why use ARROWPACK pT-Tester

- ▶ Double parametric test.
- ▶ Easy to use and very simple to read.
- ▶ For all types of rotary sealers.
- ▶ Results are instantly given by clear and simple information.
- ▶ Precise parameter measurements.
- ▶ Determines, if pressure is within optimal values.
- ▶ Determines, if temperature is within optimal values.
- ▶ Ensures integrity of closure with clear and precise interpretation.
- ▶ Combine with ARROWPACK Push-DyeTest to give complete seal assurance.



	Cat. Nr.	Qty/box
pT-Tester Steam Packing	44200	10
pT-Tester Plasma Packing	44210	10
pT-Tester Ultra	44220	10