



hm RS

Folienrollenspender für Sterilgutverpackungen

Film roll dispenser for sterile packaging

hawo

hm RS

Folienrollenspender mit Schneidevorrichtung für Sterilgutverpackungen


Der Folienspender für Sterilgutverpackungsrollen ist eine sinnvolle Ergänzung zu den hawomed- und hawodent- Foliensiegelgeräten. Die integrierte Schneidevorrichtung und die Möglichkeit, den Folienspender mit mehreren Schlauch- und Beutelfolienrollen zu bestücken, machen den hm RS zu einem unverzichtbaren Bestandteil Ihrer Pack- und Siegelstation.

Zubehör

hawomed Folienspender **hm 632 R**

Technische Daten

Ausführung	Edelstahl
hm RS BxHxT	873 x 430 x 420 mm
hm 631 S BxHxT	873 x 125 x 420 mm
hm 632 R BxHxT	858 x 425 x 200 mm
Rollendurchmesser	max. 300 mm
Schneidelänge	max. 750 mm
hm 621 S BxHxT	458 x 125 x 344 mm
Rollendurchmesser	max. 200 mm
Schneidelänge	max. 400 mm



hm 631 S + hm 632 R

Film roll dispenser with cutting device for sterile packaging


The film dispenser for sterile packaging rolls is a useful addition to the hawomed and hawodent film sealing units. The integrated cutting device and the possibility of fitting several tubular and pouch film rolls to the film dispenser make the hm RS an indispensable part of your packaging and sealing station.

Accessories

hawomed film dispenser **hm 632 R**

Technical data

Finish	Stainless steel
hm RS w x h x d	873 x 430 x 420 mm (34.3 x 16.9 x 16.5 in.)
hm 631 S w x h x d	873 x 125 x 420 mm (34.3 x 4.9 x 16.5 in.)
hm 632 R w x h x d	858 x 425 x 200 mm (33.7 x 16.7 x 7.9 in.)
Roll diameter	max. 300 mm (11.8 in.)
Cutting length	approx. 750 mm (29.5 in.)
hm 621 S w x h x d	458 x 125 x 344 mm (18.0 x 4.9 x 13.5 in.)
Roll diameter	max. 200 mm (7.9 in.)
Cutting length	max. 400 mm (15.8 in.)



hm 631 S

hawomedpack

Schlauchfolien/Tubular films

hmp 050	50 mm x 200 m (2.0 in. x 656 ft.)
hmp 075	75 mm x 200 m (3.0 in. x 656 ft.)
hmp 100	100 mm x 200 m (3.9 in. x 656 ft.)
hmp 150	150 mm x 200 m (5.9 in. x 656 ft.)
hmp 200	200 mm x 200 m (7.9 in. x 656 ft.)



Technische Änderungen vorbehalten
Subject to technical alterations

Oktober 2009 9.334.051 Printed in Germany
10/09 / 0 / WRB



hawo

hawo GmbH · Obere Au 2-4 · 74847 Obrigheim/Germany
T +49 (0)6261 9770-0 · F +49 (0)6261 9770-69
info@hawo.com · www.hawo.com