

# Rocket Craft™ Suction Pumps

The Rocket Craft™ and the Rocket DUO-VAC™ Suction Pumps have been the gold standard for reliable oocyte aspiration for over 30 year and have been developed to provide smooth, low volume vacuum at pre-determined negative pressures, so allowing simple, low turbulence, oocyte recovery.

- **HIGHLY RELIABLE:** Rocket Craft™ Pumps have a proven 30 year history of reliability and excellent value.
- **SIMPLY, SAFE OPERATION** with air-operated foot controls. High vacuum foot control available on Duo-Vac Suction Pumps.
- **NEW FOOTSWITCH:** A new style foot switch base replaces the traditional chromed plates and now include an lighter weight connection set to improve ease of handling and reliability.
- **SUPER QUIET** diaphragm pump makes it easy to live with in the OR and needs only annual service.
- **FULLY VARIABLE VACUUM** from 0-400mmHg.
- **SINGLE USE WATER TRAP SET** provides cost effective solution to non-sterile water traps and tubing.



## Ordering Information

<b>Rocket Craft™ Suction Pump</b>	240V/50Hz	<b>R29655</b>
<b>Rocket Craft™ Suction Pump</b>	120V/60Hz	<b>R29654</b>
<b>Rocket Craft™ DUO-VAC Suction Pump</b>	240V/50Hz	<b>R29660</b>
<b>Rocket Craft™ DUO-VAC Suction Pump</b>	120V/60Hz	<b>R29661</b>
<b>Rocket Water Trap Set</b> suitable for all Craft™ and GenX Suction Pumps packed individually, 10 per carton	Sterile each	<b>R57685</b>

